


MEDI-CRYL

TWO LAYER ANTERIOR AND POSTERIOR
ACRYLIC TEETH



ISO 13485:2016

CE



MR. Dental[®]
England



Upper Anterior

S13

- (1) 42.3
- (2) 7.4
- (3) 8.9
- (4) 8.4



S14

- (1) 42.6
- (2) 7.7
- (3) 9.2
- (4) 7.9




S15

- (1) 41.0
- (2) 8.6
- (3) 10.3
- (4) 7.7



S18

- (1) 45.2
- (2) 9.0
- (3) 11.0
- (4) 8.7



S25

- (1) 43.8
- (2) 11.2
- (3) 11.2
- (4) 8.2



S32

- (1) 42.8
- (2) 9.0
- (3) 10.7
- (4) 7.8



S33

- (1) 46.5
- (2) 10.0
- (3) 11.5
- (4) 8.6



S38

- (1) 41.4
- (2) 9.2
- (3) 11.1
- (4) 7.8



S39

- (1) 43.3
- (2) 9.1
- (3) 11.2
- (4) 8.0



S41

- (1) 44.4
- (2) 9.3
- (3) 11.4
- (4) 8.1



S42

- (1) 43.4
- (2) 9.0
- (3) 11.2
- (4) 8.0



S51

- (1) 39.7
- (2) 7.7
- (3) 10.0
- (4) 7.2



S56


- (1) 47.0
- (2) 10.0
- (3) 12.2
- (4) 8.1



S67

- (1) 50.4
- (2) 10.3
- (3) 12.2
- (4) 9.2






Lower Anterior

I13

- (1) 33.3
- (2) 7.4
- (3) 9.2
- (4) 5.0



I18

- (1) 36.3
- (2) 8.6
- (3) 10.4
- (4) 5.3



I38

- (1) 31.7
- (2) 8.0
- (3) 10.0
- (4) 4.7



I46

- (1) 35.1
- (2) 10.2
- (3) 10.4
- (4) 5.0




I51

- (1) 31.2
- (2) 7.3
- (3) 9.3
- (4) 5.5



I67

- (1) 40.7
- (2) 10.3
- (3) 12.7
- (4) 5.7



Upper & Lower Posterior

D35 (1) 30.7 (2) 33.0 (3) 8.5 (4) 5.5



D37 (1) 30.7 (2) 33.2 (3) 8.4 (4) 8.8



D39 (1) 30.1 (2) 35.0 (3) 8.5 (4) 8.5



D45 (1) 31.0 (2) 33.3 (3) 8.8 (4) 8.6



Articulation Combination Chart

UPPER ANTERIORS	LOWER ANTERIORS	POSTERIORS
S13	L13	D35
S14	L13	D35
S15	L38	D37
S18	L18	D39
S25	L18	D39
S32	L41	D39
S33	L46	D45
S38	L38	D37
S39	L51	D35
S41	L41	D39
S42	L41	D39
S51	L51	D35
S56	L46	D45
S67	L67	D45



Mfg by: MR. Dental
4 Manor Way, Old Woking
Surrey GU22 9JX, England
www.mrdental.co.uk



MADE IN ENGLAND